

SLC24A4 Polyclonal Antibody



Catalog Number: 166096

Product name

SLC24A4 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SLC24A4

Preparation

Antigen: A synthetic peptide of human SLC24A4

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:500 - 1:1000

Validations

Human lung



Western blot - SLC24A4 Polyclonal Antibody

Western blot analysis of extracts of human lung, using SLC24A4 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.