

# ABCD1 antibody



Catalog Number: 107096

## Product name

ABCD1 antibody

## Specificity

Human; other species not tested.

## Antibody description

ABCD1 Mouse Monoclonal antibody. Positive WB detected in Hela cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 75 kDa

## Preparation

This antibody was obtained by immunization of ABCD1 recombinant protein (Accession Number: NM\_000033). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107556(ABCD1 antibody) at dilution of 1:400



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:107556(ABCD1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)