

EWS antibody

Catalog Number: 110496

Product name

EWS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EWS Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human testis tissue, K-562 cells, mouse testis tissue. Observed molecular weight by Western-blot: 90 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005243). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110496(EWS antibody) at dilution of 1:600