

EWS antibody



Catalog Number: 110496

Product name

Polyclonal

EWS antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

EWS Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human testis tissue, K-562 cells, mouse testis tissue. Observed molecular weight by Western-blot: 90 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005243). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110496(EWS antibody) at dilution of 1:600