

Anti-Transglutaminase 2 Rabbit antibody



Catalog Number: 167496

Product name

Rabbit IgG

Anti-Transglutaminase 2 Rabbit antibody

Applications

Specificity

WB, IHC-P, IP

Human, Mouse

Dilutions

Antibody description

WB: 1/2000-1/10000

Transglutaminase 2 Rabbit polyclonal antibody

IHC: 1/50

Preparation


IP: 1/20

Antigen: A synthetic peptide of human Transglutaminase 2

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Western blot detection of Transglutaminase 2 in Rat Brain lysates using Transglutaminase 2 Rabbit pAb(1:1000 diluted). Predicted band size:77kDa.Observed band size:77kDa.

Ig Type