

COX6C antibody



Catalog Number: 109496

Product name

COX6C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

COX6C Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 9kd

Preparation

This antibody was obtained by immunization of COX6C recombinant protein (Accession Number: XM_017013020). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109496(COX6C antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:109496(COX6C antibody) at dilution of 1:50 (under 10x lens)