

# PIGR antibody

Catalog Number: 113896

## Product name

PIGR antibody

## Specificity

Human; other species not tested.

## Antibody description

PIGR Rabbit Polyclonal antibody. Positive IHC detected in human small intestine tissue, human kidney tissue.

## Preparation

This antibody was obtained by immunization of PIGR recombinant protein (Accession Number: NM\_002644). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:113896(PIGR Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:113896(PIGR Antibody) at dilution of 1:200 (under 40x lens).