

# Anti-Beclin 1 Rabbit antibody



Catalog Number: 168896

## Product name

Anti-Beclin 1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Beclin 1 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Beclin 1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

1:1000

1:20

1:20

## Validations



Western blot detection of Beclin 1 in K562,Rat Brain,C6,3T3,Hela cell lysates using Beclin 1 Rabbit mAb(1:1000 diluted).Predicted band size:52kDa.Observed band size:60kDa.