

# CBFA2T2 antibody

Catalog Number: 108896

## Product name

CBFA2T2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CBFA2T2 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in HEK-293 cells, mouse testis tissue, PC-3 cells. Observed molecular weight by Western-blot: 85kd

## Preparation

This antibody was obtained by immunization of CBFA2T2 recombinant protein (Accession Number: BC015066). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108896(CBFA2T2 antibody) at dilution of



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:108896(CBFA2T2 antibody) at dilution of 1:50 (under 10x lens)