

EIF3H antibody



Catalog Number: 110196

Product name

Polyclonal

EIF3H antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

EIF3H Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 40kd

Dilutions

Recommended Dilution:

Preparation


WB: 1:200-1:2000

This antibody was obtained by immunization of EIF3H recombinant protein (Accession Number: NM_003756). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110196(EIF3H antibody) at dilution of 1:500