

GRASP65,GORASP1 antibody



Catalog Number: 111196

Product name

GRASP65,GORASP1 antibody

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

GRASP65,GORASP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 46±65 kDa

Preparation

This antibody was obtained by immunization of GRASP65,GORASP1 recombinant protein (Accession Number: BC008928). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF, WB

Dilutions

Recommended Dilution:


WB: N/AIF: 1:10-1:100

Validations



Blocking: 5% milk in TBS-T
Primary ab: GRASP65 (1:1000) in 1%BSA in TBS-T O/N at 4 degree
Sec ab: Anti Rabbit (1:10,000) in 0.5%BSA in TBS-T for 35 min atRT
Dectction by ECL GE: 45 sec exposure:

WB result of GRASP65 anitbody (Catalog No:111196,1:1000) with HeLa cell lysate by Dr. Neeraj Tiwari, James Rothman Lab.



Immunofluorescent analysis of HeLa cells using Catalog No:111196(GRASP65,GORASP1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG