

Anti-BRG1 (6D7) Mouse antibody



Catalog Number: 168196

Product name

Anti-BRG1 (6D7) Mouse antibody

Specificity

Human, Mouse

Antibody description

BRG1 (6D7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human BRG1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF, IP

Dilutions


WB: 1/1000

;ICC: 1/50


Validations



Immunoprecipitation analysis of K562 cell lysates using BRG1 mouse mAb (168196).



Immunofluorescent analysis of HeLa cells fixed with 4% Paraformaldehyde and using anti-BRG1 mouse mAb (dilution 1:50). DAPI was used to stain



nucleus(blue).

Western blot detection of BRG1 in HeLa NE,F9 and K562 cell lysates using BRG1 mouse mAb (1:1000 diluted).Predicted band size:220KDa.Observed band size:220KDa.