

# POLR3A Polyclonal Antibody



Catalog Number: 160596

## Product name

POLR3A Polyclonal Antibody

## Specificity

Human, Rat

## Antibody description

Polyclonal antibody to POLR3A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human POLR3A (NP\_008986.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200

## Validations



Western blot - POLR3A Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using POLR3A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.