

SLIT2 Polyclonal Antibody



Catalog Number: 162596

Product name

SLIT2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SLIT2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1330-1529 of human SLIT2 (NP_004778.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - SLIT2 Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using SLIT2 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.Detection: ECL Basic Kit .Exposure time: 90s.