

PARP9 antibody



Catalog Number: 113596

Product name

PARP9 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

PARP9 Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in HeLa cells, K-562 cells, MCF7 cells, rat skeletal muscle tissue. Positive IHC detected in human spleen tissue, human brain tissue, human heart tissue, human kidney tissue, human lung tissue, human lymphoma tissue, human ovary tissue, human skin tissue, human testis tissue, human tonsillitis tissue. Observed molecular weight by Western-blot: 88kd

Preparation

This antibody was obtained by immunization of PARP9 recombinant protein (Accession Number: NM_001387875). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113596(PARP9 antibody) at dilution of 1:500




Immunohistochemical of paraffin-embedded human spleen using Catalog No:113596(PARP9 antibody) at dilution of 1:50 (under 10x lens)

PARP9 antibody

Catalog Number: 113596



Immunohistochemical of paraffin-embedded human spleen using Catalog No:113596(PARP9 antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-PARP9 (IP:Catalog No:113596, 4ug; Detection:Catalog No:113596 1:500) with MCF-7 cells lysate 4000ug.