

Anti-LXR alpha Rabbit antibody



Catalog Number: 168596

Product name

Rabbit IgG

Anti-LXR alpha Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

1:1000

LXR alpha Rabbit monoclonal antibody

1:20

Preparation

1:20

Antigen: A synthetic peptide of human LXR alpha

1:20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Validations


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of LXR alpha in Rat Brain lysates using LXR alpha Rabbit mAb(1:1000 diluted).Predicted band size:50kDa.Observed band size:50kDa.