

CSHL1 antibody



Catalog Number: 109596

Product name

CSHL1 antibody

Specificity

human; other species not tested.

Antibody description

CSHL1 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Positive IHC detected in human placenta. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of CSHL1 recombinant protein (Accession Number: NM_001318). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:109596(CSHL1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human placenta using Catalog No:109596(CSHL1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:109596(CSHL1 antibody) at dilution of 1:100 (under 40x lens)