

KIAA1602 antibody



Catalog Number: 111996

Product name

KIAA1602 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

KIAA1602 Rabbit Polyclonal antibody. Positive IHC detected in human skin. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 145-150kd

Preparation

This antibody was obtained by immunization of KIAA1602 recombinant protein (Accession Number: BC110599). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111996(KIAA1602 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human skin using Catalog No:111996(KIAA1602 antibody) at dilution of 1:50 (under 10x lens)

Immunohistochemical of paraffin-embedded human skin using Catalog No:111996(KIAA1602 antibody) at dilution of 1:50 (under 40x lens)