

MYOT antibody



Catalog Number: 112996

Product name

MYOT antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MYOT Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of MYOT recombinant protein (Accession Number: NM_006790). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112996(MYOT antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:112996(MYOT antibody) at dilution of 1:100 (under 10x lens)