

Anti-GPR137C antibody



Catalog Number: 175996

Product name

Anti-GPR137C antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to GPR137C

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gpr137c. aa range:256-305

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.;; Western Blot analysis of various cells using GPR137C Polyclonal Antibody