

# Anti-Thrombomodulin Rabbit antibody



Catalog Number: 168996

## Product name

Anti-Thrombomodulin Rabbit antibody

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Specificity

Human

## Dilutions

1:2000

## Antibody description

Thrombomodulin Rabbit monoclonal antibody

1:50

1:50

## Preparation

Antigen: A synthetic peptide of human Thrombomodulin


1:20

1:20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations



## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

Western blot detection of Thrombomodulin in Jurkat cell lysates using Thrombomodulin Rabbit mAb(1:1000 diluted).Predicted band size:60kDa.Observed band size:95kDa.