

IGHM antibody

Catalog Number: 111696

Product name

IGHM antibody

Specificity

Human; other species not tested.

Antibody description

IGHM Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 75 kDa

Preparation

This antibody was obtained by immunization of IGHM recombinant protein (Accession Number: BC009851). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:111696(IGHM Antibody) at dilution of 1:300