

Mouse S100A11 / S100C (His Tag) recombinant protein



Catalog Number: 502696

General Information

Gene Name Synonym

Calgizzarin; Endothelial monocyte-activating polypeptide; Protein S100-C; S100 calcium-binding protein A11

Protein Construction

A DNA sequence encoding the mouse S100A11 (NP_058020.1) (Met1-Ile98) was expressed with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant mouse S100A11 consists of 114 amino acids and predicts a molecular mass of 13.2 KDa. It migrates as an approximately 14 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, 10% glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse S100A11 / S100C Protein (His Tag) SDS-PAGE