

# Anti-TFPI Rabbit antibody



Catalog Number: 167696

## Product name

Anti-TFPI Rabbit antibody

## Specificity

Human

## Antibody description

TFPI Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human TFPI

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC, IP


## Dilutions

WB: 1/1000

FC: 1/50

IP: 1/20

## Validations



Western blot detection of TFPI in A549,HL-60,U2OS cell lysates using TFPI Rabbit mAb(1:1000 diluted).Predicted band size:35kDa.Observed band size:40-55kDa.