

# Anti-FOXO4 Rabbit antibody



Catalog Number: 168696

## Product name

Anti-FOXO4 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FOXO4 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human FOXO4

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC, IP

## Dilutions


1:1000-1:5000

1:20-1:50

1:20-1:100

1:20-1:50

## Validations



Western blot detection of FOXO4 in K562,C6,3T3 cell lysates using FOXO4 Rabbit mAb(1:1000 diluted).Predicted band size:54kDa.Observed band size:65kDa.