

Human CRNN/C1orf10/DRC1/PDRC1/SEP53 (His tag) recombinant protein



Catalog Number: 518696

General Information

Gene Name Synonym

Cornulin, 53 kDa Putative Calcium-Binding Protein, 53 kDa Squamous Epithelial-Induced Stress Protein, 58 kDa Heat Shock Protein, Squamous Epithelial Heat Shock Protein 53, Tumor-Related Protein, CRNN, C1orf10, DRC1, PDRC1, SEP53

Protein Construction

A DNA sequence encoding the Human CRNN/C1orf10/DRC1/PDRC1/SEP53 Met1-Ser140 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human CRNN/C1orf10/DRC1/PDRC1/SEP53 (His tag) recombinant protein