

TYR Polyclonal Antibody



Catalog Number: 163396

Product name

TYR Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TYR

Preparation

Antigen: Recombinant protein of human TYR.

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF, IP

Dilutions


WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

IP 1:50 - 1:200

Validations



Immunofluorescence - TYR Polyclonal Antibody

Immunofluorescence analysis of rat skin using TYR antibody at dilution of 1:100. Blue: DAPI for nuclear staining.