

Anti-ZNHIT1 antibody



Catalog Number: 176396

Product name

Anti-ZNHIT1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to ZNHIT1

Preparation

This antigen of this antibody was synthesized peptide derived from znhit1 . at aa range: 40-120

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of various cells using ZNHIT1 Polyclonal Antibody diluted at 1: 1000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000