

# Anti-ERM/ETV5 (3H4) Mouse antibody



Catalog Number: 168396

## Product name

Anti-ERM/ETV5 (3H4) Mouse antibody

WB: 1/500 - 1/2000

IF: 1/200- 1/500

## Specificity

Human, Rat, Mouse

ELISA: 1/10000

## Antibody description

ERM/ETV5 (3H4) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human ETV5 expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type


Mouse IgG1

## Applications


WB, ICC/IF

## Dilutions

## Validations



Immunocytochemistry staining of HeLa cells fixed with 1% Paraformaldehyde and using anti-ERM/Etv5 mouse mAb (dilution 1:500).



Western blot detection of ERM/Etv5 in C6,3T3 and K562 cell lysates using ERM/Etv5 mouse mAb (1:500 diluted). Predicted band size:58KDa. Observed band size:58KDa.