

# PGT antibody



Catalog Number: 113796

---

## Product name

Polyclonal

PGT antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

PGT Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, mouse heart tissue, mouse liver tissue. Observed molecular weight by Western-blot: 65 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of PGT recombinant protein (Accession Number: NM\_005630). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113796(PGT antibody) at dilution of 1:1000

## Clonality