

# RPAP1 antibody



Catalog Number: 114796

## Product name

RPAP1 antibody

## Immunogen

[Human RPAP1 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse; other species not tested.

## Antibody description

RPAP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse brain tissue.

## Preparation

This antibody was obtained by immunization of RPAP1 recombinant protein (Accession Number: NM\_015540). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114796(RPAP1 antibody) at dilution of 1:800