

# Anti-ACE2 Rabbit antibody



Catalog Number: 166796

## Product name

Anti-ACE2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

ACE2 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human ACE2

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


WB: 1/1000-1/2000

IHC: 1/50-1/200

ICC/IF: 1/50-1/200


IP: 1/50

## Validations



Immunohistochemical analysis of paraffin-

embedded human kidney



Western blot analysis of ACE2 expression in Human kidney lysate.