

# ADNP antibody



Catalog Number: 107796

## Product name

ADNP antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ADNP Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 135-145 kDa

## Preparation

This antibody was obtained by immunization of ADNP recombinant protein (Accession Number: NM\_001347511). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107796(ADNP antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells, using ADNP antibody Catalog No:107796 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).